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*Universal adaptor
for Arrow pump*



*Universal adaptor
for DataScope pump*


Datascope	Arrow
97,98, 98XT, CS100, CS300	ACAT, AutoCat, AutoCat 2, AutoCat 2 Wave, KAAT II

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
8Fr Response Performance in a 7Fr device

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

8Fr Response Performance in a 7Fr device⁽²⁾

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
The fiber-optic sensor instantaneously transmits a precise pressure signal to the CS300

- Reduced signal transfer delay vs. fluid-filled column
- Maximizes helium gas shuttle lumen
- Speeds and standardizes pressure set-up
- **Just plug it in and go!**
- Advanced tip design and reduced signal artifact
- Provides a clean, crisp arterial pressure waveform
- Arterial pressure waveform not affected by electro-surgical interference or patient movement and positioning

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doi:10.1510/icvts.2009.205823

INTERACTIVE
CARDIOVASCULAR AND
THORACIC SURGERY

Interactive CardioVascular and Thoracic Surgery 9 (2009) 135-137

www.icvts.org

Case report - Assisted circulation

Safety and efficacy of transbrachial intra-aortic balloon pumping with the use of 7-Fr catheters in patients undergoing coronary bypass surgery[☆]

Antonino S. Rubino*, Francesco Onorati, Filiberto Serraino, Attilio Renzulli

Cardiac Surgery Unit, Magna Graecia University of Catanzaro, Viale Europa, Località Germaneto, 88100 Catanzaro, Italy

Received 21 February 2009; received in revised form 1 April 2009; accepted 2 April 2009

Abstract

We report the cases of five consecutive patients undergoing coronary artery bypass grafting (CABG) who required a transbrachial approach for 7-Fr catheter intra-aortic balloon pumping (IABP) insertion because of unsuitable femoral arteries. No adverse outcomes occurred in any patient during a mean 72 h of IABP support. Our experience with 7-Fr catheters appears to confirm previous reports of the safety and efficacy of transbrachial IABP assistance and suggests that such support can be provided safely for an extended duration with the use of these smaller catheters.

© 2009 Published by European Association for Cardio-Thoracic Surgery. All rights reserved.

Keywords: Coronary artery bypass grafting; Medical technology; Assisted circulation; Counterpulsation; Intra-aortic balloon pumping

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doi:10.1510/icvts.2009.207522

INTERACTIVE
CARDIOVASCULAR AND
THORACIC SURGERY

Interactive CardioVascular and Thoracic Surgery 9 (2009) 369-371

www.icvts.org

Case report - Assisted circulation

Intra aortic balloon pump insertion through left axillary artery in patients with severe peripheral arterial disease

Giuseppe Zattera^a, Pasquale Totaro^{a,*}, Andrea Maria D'Armini^{a,b}, Mario Viganò^{a,b}

^a*Division of Cardiac Surgery, IRCCS Foundation San Matteo, Pavia, Italy*

^b*School of Medicine, University of Pavia, Italy*

Received 17 March 2009; received in revised form 15 April 2009; accepted 20 April 2009

Abstract

Intra aortic balloon pump (IABP) is the mechanical assist device most frequently used in cardiac surgery. Recent demonstration of better outcome following preoperative IABP insertion in high-risk patients has further extended its indication. However, due to an increasing complexity of patients currently referred for cardiac surgery, several patients with potential indication for preoperative and/or postoperative IABP present severe peripheral vascular disease which usually contraindicates IABP insertion. Here we present an alternative technique for IABP insertion in patients with severe peripheral vessel disease.

© 2009 Published by European Association for Cardio-Thoracic Surgery. All rights reserved.

Keywords: Cardiac surgery; Intra aortic balloon pump; Vascular complications

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Subclavian Artery Access for Ambulatory Balloon Pump Insertion

Jai Raman, MD, PhD, Gabriel Loor, MD, Mark London, MD, and Neeraj Jolly, MD, DM

Sections of Cardiothoracic Surgery and Cardiology, University of Chicago Medical Center, Chicago, Illinois

Intraaortic balloon pump insertion is traditionally performed through the femoral artery in the groin. However, this restricts the patient to bed rest, and prolonged implantation can be associated with infections in the groin crease. We describe a technique of insertion of a balloon pump through the subclavian artery, which al-

lows the patient to ambulate. This technique can also be performed under local anesthesia in the cardiac catheterization laboratory.


(Ann Thorac Surg 2010;90:1032–4)

© 2010 by The Society of Thoracic Surgeons

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ULTRA - IABP Patient Sizing Chart

7Fr Intra Aortic Balloon Catheter

See instructions for use for full list of compatible pumps.

Federal Law (USA) restricts this device to sale by, or on the order of, a physician.

INSIGHTRA
MEDICAL

9200 Irvine Center Drive, Ste. 200
Irvine, CA 92618 USA • www.insightra.com

© 2010 Insightra Medical Inc. USA. All trademarks are owned by their respective companies

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
- **Comparative Study of Five Types of IABP balloons in terms of incidence of balloon rupture and other complications: A multi institutional study**
Nishida et al : Artificial organs 18(10) 746-751 1994
- 1988 - 1992 (9Fr with equivalent components to 7Fr)
- 2,876 patients, 14 hospitals in Japan
- Looking at catheter related complications
- Results

Manufacturer					
Complication	TMP	KONTRON	DATASCOPE	BOST.SCI	ST.JUDE
Balloon rupture	0	15	20	3	11
Other	26	41	9	6	7
Total %	26(2.7%)	56(5.9%)	29(6.0%)	9 (4.0%)	18 (9.6%)

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- **Clinical use evaluation of the TMP 7-Fr IABP balloon catheter**

Matsuda et al Jan 2001 - May 2005

- 280 patients to compare TMP 7Fr to (Datascope, Zeon, 8Fr TMP and 10Fr TMP and 10Fr TMP)
First published - Cardioangiology 2002;51:465–466. 23
- Evaluated: Percutaneous insertion and maintain position, antiarrhythmia responsiveness, pumping performance, anti-ischemic effects and compatibility with Datascope and Arrow IAB pumps.
- Results:
 - The TMP 7Fr balloon had the lowest insertion resistance
 - All products stayed in position with maximal change of 4mm
 - Patients with Ventricular tachicardia(VT) and Afib were examined for antiarrhythmia responsiveness. All products comparable
 - Anti-ischemia. No reports of lower limb ischemia with the 7Fr
 - Two cases of thrombus formation - Not with the TMP products
 - 7Fr worked well with both Datascope & Arrow pumps and compared to the Datascope 8Fr Catheter performance (Plateau pressures and inflation/deflation time were in-line with the 8Fr catheters
 - Pressure monitoring was equally as good
- **Summary: The 7Fr balloon performed as well as the 8Fr balloons (all manufacturers) in terms of performance, but were easier to insert.**

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Optimization and safety of the Intra-Aortic balloon pumping catheter

Tsutsui et al

J Artif organs (2007) 10,143:148

Summary article talking of the background development of the TMP devices

The range of devices available

Details are given of the way they managed to reduce diameter size without losing performance

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Effects of Intra-Aortic Balloon Pumping on Graft Flow in Coronary Surgery: An Intraoperative Transit-Time Flowmetric Study

Yoshiyuki Takami, MD, and Hiroshi Masumoto, MD

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Prolonged Use for at Least 10 Days of Intraaortic Balloon Pumping (IABP) for Heart Failure

Kiyohiro OSHIMA,^{1,2} MD, Yasuo MORISHITA,¹ MD, Hiroshi HINOHARA,² MD, Yuji KADOI,² MD, Yoshiro HAYASHI,² MD, Yukio TAJIMA,³ ME, and Fumio KUNIMOTO,² MD

Intl Heart J Vol 46 No 6 November 2005

18 patients underwent extended (10 day) therapy using the 7Fr IABP.

No reports of balloon rupture even at extended use

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