


## EZ Glide Aortic Perfusion Cannula

The EZ Glide aortic perfusion cannula is intended to create a dispersive flow.

- Tip design disperses return flow in a conical spray pattern
- Unique auto-dilating tip

**10 units per case**



## Straight Cannula

**14 inch (35 cm) overall length**

3/8 inch vented connector

EZS21TA 21 Fr (7.0 mm)

EZS24TA 24 Fr (8.0 mm)



3/8 inch non-vented connector

EZS21A 21 Fr (7.0 mm)

EZS24A 24 Fr (8.0 mm)



## EZ Glide Aortic Perfusion Cannula (*continued*)

### Curved Cannula with Suture Bump

---

14.8 inch (37.6 cm) overall length

3/8 inch vented connector

EZC21TA 21 Fr (7.0 mm)

EZC24TA 24 Fr (8.0 mm)



3/8 inch non-vented connector

EZC21A 21 Fr (7.0 mm)

EZC24A 24 Fr (8.0 mm)



### Curved Cannula with Suture Flange

---

14.8 inch (37.6 cm) overall length

3/8 inch vented connector

EZF21TA 21 Fr (7.0 mm)


EZF24TA 24 Fr (8.0 mm)



3/8 inch non-vented connector

EZF21A 21 Fr (7.0 mm)

EZF24A 24 Fr (8.0 mm)



## Arterial Cannula Accessories

These accessories can be used in conjunction with Edwards arterial cannulae.

### Vascular Tourniquet

2 per pouch

---

TK2 5 inch (12.7 cm) sheath

