

NO	Namw	QTY	Drawing	Package	Remark
1		1		Blue Wrap + Breath Bag	
2		1			
3		1			
4	3				
5		3			
6		1		Out packing: paper breathing bag, paper tray and breath bag.	
7		1			
8		1			
9	3				
10		1			
11		1			
12		1			
13		1			
14		1			
15		1			
16	上	10		Individual package for connector 150x175 peel bag	
17		1			
18		1			
19		1			
20		5	5PCS WIDE TIE BANDS	Individual package 75*260 peel bag.	
21		1	20 ml SYRINGE WITH LUER PORT	Individual package 150*175 peel bag.	

1. Wrap 1-3 with blue wrap (900*900mm) and then put it into breathing bag in 370*560mm.
 2. coil every line of line 6-12.
 3. put pack 1 and 2 into paper tray. Number 1 is on top and pack 2 is at the bottom. Then pack it by 1200*1200biowrap. And paste label and E0 sticker.
 4. put pack 3 into 530*690 breath bag. Then seal it and paste carton label.

